Readiness for Diversity

Pressing Schedules

Horizon (Shenpaz)¹

Translucency	Size	Shade	Investment Ring	Start Temperature	Heating Rate	Max Temperature	Holding Time	Vacuum On	Vacuum Off
НТ		A1, A2, A3, A3.5, B1, B2, B3, B4, C1, C2, C3, C4, D2, D3, D4, W1, W2, W3, W4		g) ge 700°C 60	60℃/min	915℃	15 Min / 20 Min	700℃	915℃
LT	R10 / R20	' ' ' '	Small (100g) Large (200g)						
МО		M00, M01, M02, M03, M04				920℃			920℃

Austromat Press-i-dent (Dekema)²

Translucency	Size	Shade	Investment Ring	Start Temperature	Heating Rate	Final Temperature	Holding Time	Press Duration	Press Level
НТ		A1, A2, A3, A3.5, B1, B2, B3, B4, C1, C2, C3, C4, D2, D3, D4, W1, W2, W3, W4				925℃			
LT	R10 / R20	A1, A2, A3, A3.5, B1, B2, B3, B4, C1, C2, C3, C4, D2, D3, D4 W1, W2, W3, W4	Small (100g) / Large (200g)	700℃	60°C/min	9230	20 Min (100g) / 30 Min (200g)	Auto 1	6
МО		M00, M01, M02, M03, M04				930℃			

EP600 (Ivoclar Vivadent)³

Translucency	Size	Shade	Investment Ring	Stand-by Temperature	Temperature Increase	Holding Temperature	Holding Time	Stop Speed
НТ		A1, A2, A3, A3.5, B1, B2, B3, B4, C1, C2, C3, C4, D2, D3, D4, W1, W2, W3, W4						
LT	R10 / R20	A1, A2, A3, A3.5, B1, B2, B3, B4, C1, C2, C3, C4, D2, D3, D4 W1, W2, W3, W4	Small (100g) / Large (200g)	700℃	60°C/min	930℃	15 Min (100g) / 25 Min (200g)	300 <i>µ</i> m/min
МО		M00, M01, M02, M03, M04						

Pressing Schedules

EP3000 (Ivoclar Vivadent)³

Translucency	Size	Shade	Investment Ring	Stand-by Temperature	Temperature Increase	Holding Temperature	Holding Time	Stop Speed
НТ		A1, A2, A3, A3.5, B1, B2, B3, B4, C1, C2, C3, C4, D2, D3, D4, W1, W2, W3, W4						
LT	R10 / R20	A1, A2, A3, A3.5, B1, B2, B3, B4, C1, C2, C3, C4, D2, D3, D4 W1, W2, W3, W4	Small (100g) / Large (200g)	700℃	60°C/min	915℃	15 Min (100g) / 25 Min (200g)	300µm/min
МО		M00, M01, M02, M03, M04						

EP5000 (Ivoclar Vivadent)3

Translucency	Size	Shade	Investment Ring	Stand-by Temperature	Temperature Increase	Holding Temperature	Holding Time	Stop Speed
НТ		A1, A2, A3, A3.5, B1, B2, B3, B4, C1, C2, C3, C4, D2, D3, D4, W1, W2, W3, W4						
LT	R10 / R20	A1, A2, A3, A3.5, B1, B2, B3, B4, C1, C2, C3, C4, D2, D3, D4 W1, W2, W3, W4	Small (100g) / Large (200g)	700℃	60°C/min	915℃	20 Min (100g) / 30 Min (200g)	300 <i>µ</i> m/min
МО		M00, M01, M02, M03, M04						

NOTE: The above schedules are referential guideline only

There may be a difference between the displayed temperature and the real temperature of each furnace.

When you use the Amber ingots, please verify the above standard schedule is suitable for your press furnace.

If it is not, please try to find the optimum temperature through the following process.

- 1) If there are some traces of tiny bubble on the surface of the restoration
 - \Rightarrow Please reduce the maximum temperature by 5~10 $^{\circ}$ C or holding time and try pressing again.
- 2) If the marginal area of the restoration is not formed completely
 - \Rightarrow Please increase the maximum temperature by 5~10 $^{\circ}$ C or holding time and try pressing again.
- 1. Horizon is a registered trademark of Shenpaz Dental Ltd.
- 2. Austromat Press-i-dent is a registered trademark of DEKEMA Dental-Keramiköfen GmbH.
- 3. EP600, EP3000, EP5000 are registered trademarks of Ivoclar Vivadent.

HASS Corporation

77-14, Gwahakdanji-ro, Gangneung-si, Gangwon-do, KOREA 25452

Tel: +82-70-7712-1300 / Fax: +82-33-644-1231 Customer Support: +82-2-2083-1367

E-mail: hasscorp@hassbio.com Website: www.hassbio.com This material is designed for usage in dentistry. Follow instructions HASS is not liable for any loss caused by failure to comply with regulation or scope of indication. Users are responsible for testing products to verify the compatibility for any usage which are not written in the instructions. The explanations and data contained within do not carry any guarantees and/or obligations. All enclosed recommendations and restrictions apply

All enclosed recommendations and restrictions apply when used with products from other manufacturers.



Lithium Disilicate-Based Press Ingots Amber Press








Whatever You Imagine

Superior Strength

Amber® Press's flexural strength is higher than conventional Lithium Disilicate press materials. Amber Press presents the biaxial flexural strength of 460 MPa.



Simple and Safe

After pressing, very little reaction layer remains on Amber® Press.
There is no need to apply any acid for clean-up, thereby ensuring a simple and nonhazardous process.



Aesthetic Results with Amber®Press

Amber® Press helps to create highly aesthetic and natural-looking dental restorations, which provides patients beautiful and natural smile.




Restoration of 12 through 22 using Amber Press M01 by layering technique Courtesy of CDT. Seungsub Lee, Seoul, Korea

Outstanding Versatility

Compatible with Various Veneering Materials

Amber® Press ingots are compatible with various veneering materials for lithium disilicate.



- *, ** Not a registered trademark of HASS Corp.
- ** If you use Initial Zr-FS (GC) or Creation ZI-F (Creation Willi Geller), we recommend following the Amber® Press Guidelines for Zirconia veneering materials.

Indications



Inlays



Onlays



Veneers



Anterior Single Crowns



Posterior Single Crowns



3-Unit Bridge
* up to the second Premolar

Readiness for Diversity


Product Line-up


Ambe	r ® Press	Dimensions (mm)	pcs / Pack	
	R10		5 Ingots	
	R20	Ø12.7×T20	3 Ingots	

Diverse Shade Option



Amber® Press provides three levels of translucency and over 43 shades overall. It will make any restoration come to life.





MO(Medium Opacity)

